


Collagen VIII alpha 2 Rabbit pAb

ABZ1654

Product Overview

Name	Collagen VIII alpha 2 Rabbit pAb
Accession(Primary)	P25067
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P25067)

Gene Name COL8A2**About Gene Symbol**

Western blot analysis of Collagen VIII alpha 2 in hela, 3T3 lysates using Collagen VIII alpha 2 antibody.

