


GCM2 Rabbit pAb

ABZ1655

Product Overview

Name	GCM2 Rabbit pAb
Accession(Primary)	O75603
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (O75603)

Gene Name GCM2**About Gene Symbol**

Western blot analysis of GCM2 in 293T lysates
using GCM2 antibody.

