

MAP4K2/3 Rabbit pAb

ABZ1656


Product Overview

Name	MAP4K2/3 Rabbit pAb
Accession(Primary)	Q12851
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q12851)

Gene Name MAP4K2/MAP4K3

About Gene Symbol



Western blot analysis of MAP4K2/3 in 293T, hela lysates using MAP4K2/3 antibody.

