


C8 alpha Rabbit pAb

ABZ1657

Product Overview

Name	C8 alpha Rabbit pAb
Accession(Primary)	P07357
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 65 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P07357)

Gene Name C8A**About Gene Symbol**

Western blot analysis of C8 alpha in mouse kidney, mouse heart lysates using C8 alpha antibody.

