


Collagen XIII alpha 1 Rabbit pAb

ABZ1658

Product Overview

Name	Collagen XIII alpha 1 Rabbit pAb
Accession(Primary)	Q5TAT6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q5TAT6)

Gene Name COL13A1**About Gene Symbol**

Western blot analysis of Collagen XIII alpha 1 in HeLa, mouse kidney lysates using Collagen XIII alpha 1 antibody.

