

LATS1 Rabbit pAb

ABZ1659


Product Overview

Name	LATS1 Rabbit pAb
Accession(Primary)	O95835
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 127 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (O95835)

Gene Name LATS1

About Gene Symbol



1 HEPG2 UV
2 CAC02

Western blot analysis of LATS1 in various lysates using LATS1 antibody.

