

TMEM72 Rabbit pAb

ABZ1660


Product Overview

Name	TMEM72 Rabbit pAb
Accession(Primary)	A0PK05
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (A0PK05)

Gene Name TMEM72

About Gene Symbol



1 mouse-liver
2 CAC02

Western blot analysis of Transmembrane Protein
72 in various lysates using Transmembrane Protein
72 antibody.