

B3GALT2 Rabbit pAb

ABZ1661


Product Overview

Name	B3GALT2 Rabbit pAb
Accession(Primary)	O43825
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:100-1:300; IF-1:200-1:1000

Target information (O43825)

Gene Name B3GALT2

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human heart using B3GALT2 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.