

beta 4 GalNAcT1 Rabbit pAb


ABZ1662

Product Overview

Name	beta 4 GalNAcT1 Rabbit pAb
Accession(Primary)	Q8TDX6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q8TDX6)

Gene Name	CSGALNACT1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using beta 4 GalNAcT1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.