

Human kappa Light Chain Rabbit pAb

ABZ1663


Product Overview

Name	Human kappa Light Chain Rabbit pAb
Accession(Primary)	P01834
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000


Target information (P01834)

Gene Name IGKC

About Gene Symbol



Indirect ELISA assay for Rabbit Human kappa light chain antibody. Antigen coating concentration: 2ug/ml.



Western blot analysis of Human kappa Light Chain in total Human IgG reduced and denatured in 2X SDS loading buffer with DTT using Human kappa Light Chain antibody