

Human lambda Light Chain Rabbit pAb

ABZ1664


Product Overview

Name	Human lambda Light Chain Rabbit pAb
Accession(Primary)	P0CG04
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 11 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000


Target information (P0CG04)

Gene Name IGLC1

About Gene Symbol



Indirect ELISA assay for Rabbit Human kappa light chain antibody. Antigen coating concentration: 2ug/ml.



Western blot analysis of Human lambda Light Chain in total Human IgG reduced and denatured in 2X SDS loading buffer with DTT using Human lambda Light Chain antibody