

BTN3A Rabbit pAb

ABZ1665


Product Overview

Name	BTN3A Rabbit pAb
Accession(Primary)	O00481
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (O00481)

Gene Name BTN3A1/2/3

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human kidney using BTN3A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using BTN3A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.