

# Myoblast Determination Protein 1 Rabbit pAb

ABZ1669

## Product Overview

|                             |                                             |
|-----------------------------|---------------------------------------------|
| <b>Name</b>                 | Myoblast Determination Protein 1 Rabbit pAb |
| <b>Accession(Primary)</b>   | P15172                                      |
| <b>Host Species</b>         | Rabbit                                      |
| <b>Clonality</b>            | Polyclonal                                  |
| <b>Isotype</b>              | IgG                                         |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Chicken                  |
| <b>Calculated MW</b>        | Calculated MW: 35 kDa; Observed MW: 35 kDa  |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:50-1:300             |

## Target information (P15172)

**Gene Name** MYOD1

**About Gene Symbol**



