

Cytochrome B Rabbit pAb

ABZ1671


Product Overview

Name	Cytochrome B Rabbit pAb
Accession(Primary)	P00156
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 43 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P00156)

Gene Name MT-CYB

About Gene Symbol



Western blot analysis of Cytochrome B in mouse lung, mouse brain 3, mouse spleen 4, mouse kidney 5, mouse heart lysates using Cytochrome B antibody.

