

Granzyme B Rabbit pAb


ABZ1672

Product Overview

Name	Granzyme B Rabbit pAb
Accession(Primary)	P10144
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (P10144)

Gene Name	GZMB
About Gene Symbol	



Western blot analysis of Granzyme B in mouse lung lysates using Granzyme B antibody.