

CD44 Rabbit pAb


ABZ1673

Product Overview

Name	CD44 Rabbit pAb
Accession(Primary)	P16070
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:200; IHC-1:50-1:100

Target information (P16070)

Gene Name	CD44
About Gene Symbol	CD44 (CD44 Molecule (Indian Blood Group)) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify t



Western blot analysis of CD44 in various lysates using CD44 antibody.