

GRP78 BiP Rabbit pAb


ABZ1676

Product Overview

Name	GRP78 BiP Rabbit pAb
Accession(Primary)	P11021
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-P:1:50-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (P11021)

Gene Name	HSPA5
About Gene Symbol	HSPA5 (Heat Shock Protein Family A (Hsp70) Member 5) is a molecular chaperone involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mecha



Western blot analysis of GRP78 BiP in various lysates using GRP78 BiP antibody.