

SNAIL Rabbit pAb


ABZ1679

Product Overview

Name	SNAIL Rabbit pAb
Accession(Primary)	O95863
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (O95863)

Gene Name	SNAI1
About Gene Symbol	



Western blot analysis of SNAI1 in 3T3 lysates
using SNAI1 antibody.