

c-Jun Rabbit pAb


ABZ1680

Product Overview


Name	c-Jun Rabbit pAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P05412)


Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of c-Jun in Jurkat and HeLa lysates using c-Jun antibody.



Western blot analysis of c-Jun in KB lysates using AP1 antibody.



Western blot analysis of c-Jun in KB MCF-7 lysates
using AP1 antibody