

RPL7 Rabbit pAb


ABZ1682

Product Overview


Name	RPL7 Rabbit pAb
Accession(Primary)	P18124
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P18124)


Gene Name	RPL7
About Gene Symbol	




Western blot analysis of RPL7 in various lysates using RPL7 antibody.




Western blot analysis of RPL7 in cat sample lysates using Ribosomal Protein L7 antibody.



Western blot analysis of RPL7 in mouse kidney, mouse brain, HELA, KB, SH-SY5Y, 293T and 3T3 lysates using Ribosomal Protein L7 antibody.



Western blot analysis of RPL7 in 293 and HeLa lysates using RPL7 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of PL7 in K562 using RPL7 antibody.