

VEGFA Rabbit pAb

ABZ1685


Product Overview

Name	VEGFA Rabbit pAb
Accession(Primary)	P15692
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Rabbit, Pig
Calculated MW	Calculated MW: 27 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300.; IF-1:50-1:200


Target information (P15692)

Gene Name VEGFA


About Gene Symbol



Western blot analysis of VEGFA in 22RV SW480, BT474, K562, mouse liver lysates using VEGFA antibody.



Immunohistochemistry analysis of paraffin-embedded Human heart using VEGFA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human heart using VEGFA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.