

CD63 Rabbit pAb


ABZ1686

Product Overview

Name	CD63 Rabbit pAb
Accession(Primary)	P08962
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P08962)

Gene Name	CD63
About Gene Symbol	



Western blot analysis of CD63 in Jurkat lysates using CD63 antibody.