

Cyclin E1 Rabbit pAb


ABZ1687

Product Overview

Name	Cyclin E1 Rabbit pAb
Accession(Primary)	P24864
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P24864)

Gene Name	CCNE1
About Gene Symbol	CCNE1 (Cyclin E1) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin E1 in 3T3 and HeLa lysates using Cyclin E1 antibody.