

E Cadherin Rabbit pAb


ABZ1688

Product Overview


Name	E Cadherin Rabbit pAb
Accession(Primary)	P12830
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 97 kDa; Observed MW: 125-130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P12830)

Gene Name	CDH1
About Gene Symbol	



Western blot analysis of E-Cadherin in Jurkat, 3T3, COS7 and C6 lysates using E-Cadherin antibody.



Western blot analysis of E-Cadherin in T47D lysates using E-Cadherin antibody.