

## EIF2S1 Rabbit pAb


ABZ1689

### Product Overview

|                             |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| <b>Name</b>                 | EIF2S1 Rabbit pAb                                                 |
| <b>Accession(Primary)</b>   | P05198                                                            |
| <b>Host Species</b>         | Rabbit                                                            |
| <b>Clonality</b>            | Polyclonal                                                        |
| <b>Isotype</b>              | IgG                                                               |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Cat                                            |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa                        |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:10000 |

### Target information (P05198)

|                          |                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | EIF2S1                                                                                                                                                                                                                                                                                                                                     |
| <b>About Gene Symbol</b> | EIF2S1 (Eukaryotic Translation Initiation Factor 2 Subunit Alpha) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of eIF2A in various lysates using eIF2A antibody.