

Osteopontin Rabbit pAb


ABZ1690

Product Overview

Name	Osteopontin Rabbit pAb
Accession(Primary)	P10451
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:500; IF/ICC1:100-1:500; ELISA-1:5000-1:20000

Target information (P10451)

Gene Name	SPP1
About Gene Symbol	



Western blot analysis of Osteopontin in various
lysates using Osteopontin antibody.