

Claudin 5 Rabbit pAb


ABZ1692

Product Overview


Name	Claudin 5 Rabbit pAb
Accession(Primary)	O00501
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (O00501)

Gene Name	CLDN5
About Gene Symbol	



Western blot analysis of Claudin 5 in A549 lysates using Claudin 5 antibody.



Western blot analysis of Claudin 5 in HELA lysates using Claudin5 antibody.