

COL1A2 Rabbit pAb


ABZ1693

Product Overview

Name	COL1A2 Rabbit pAb
Accession(Primary)	P08123
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 129 kDa; Observed MW: 200 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P08123)

Gene Name	COL1A2
About Gene Symbol	



Western blot analysis of Collagen I alpha 2 in 823, 22RV293UV lysates using Collagen I alpha 2 antibody.