

IL-6 Rabbit pAb


ABZ1694

Product Overview

Name	IL-6 Rabbit pAb
Accession(Primary)	P05231
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P05231)


Gene Name	IL6
About Gene Symbol	




Western blot analysis of IL-6 in T47D, COS7 and HeLa lysates using IL-6 antibody.



Western blot analysis of IL-6 in K562 lysates using IL-6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using IL-6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lung using IL-6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.