

RANKL Rabbit pAb

ABZ1697

Product Overview

Name	RANKL Rabbit pAb
Accession(Primary)	O14788
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (O14788)

Gene Name	TNFSF11
About Gene Symbol	