

Neuregulin-1/SMDF Rabbit pAb

ABZ1698

Product Overview

Name	Neuregulin-1/SMDF Rabbit pAb
Accession(Primary)	Q02297-10
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (Q02297-10)

Gene Name	NRG1
About Gene Symbol	