

## CD14 Rabbit pAb

ABZ1699

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CD14 Rabbit pAb                            |
| <b>Accession(Primary)</b>   | P08571                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 50 kDa |
| <b>Recommended Dilution</b> | IHC-P-1:200-1:500; FC-1:50-1:100           |

### Target information (P08571)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | CD14 |
| <b>About Gene Symbol</b> |      |