

MARCH3 Rabbit pAb

ABZ1700

Product Overview

Name	MARCH3 Rabbit pAb
Accession(Primary)	Q86UD3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW:29kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q86UD3)

Gene Name	MARCHF3
About Gene Symbol	