

Caspase 4 Rabbit pAb

ABZ1701

Product Overview

Name	Caspase 4 Rabbit pAb
Accession(Primary)	P49662
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (P49662)

Gene Name	CASP4
About Gene Symbol	