

TSG6 Rabbit pAb

ABZ1703

Product Overview

Name	TSG6 Rabbit pAb
Accession(Primary)	P98066
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW31kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (P98066)

Gene Name	TNFAIP6
About Gene Symbol	