

## SLC38A2 Rabbit pAb

ABZ1705

### Product Overview

|                             |                                 |
|-----------------------------|---------------------------------|
| <b>Name</b>                 | SLC38A2 Rabbit pAb              |
| <b>Accession(Primary)</b>   | Q96QD8                          |
| <b>Host Species</b>         | Rabbit                          |
| <b>Clonality</b>            | Polyclonal                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat               |
| <b>Calculated MW</b>        | Calculated MW:56kDa             |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:50-1:300 |

### Target information (Q96QD8)

|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | SLC38A2 |
| <b>About Gene Symbol</b> |         |