

Cellubrevin Rabbit mAb


ABZ1706

Product Overview


Name	Cellubrevin Rabbit mAb
Accession(Primary)	Q15836
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q15836)

Gene Name	VAMP3
About Gene Symbol	



Western blot analysis of ubrevin in HeLa, A549, MCF-7, 293T lysates using ubrevin antibody.



Immunocytochemistry analysis of Cellubrevin (green) in hela using Cellubrevin antibody, and DAPI(blue)