

Phospholipase C beta 1/PLCB1 Rabbit mAb

ABZ1707

Product Overview

Name	Phospholipase C beta 1/PLCB1 Rabbit mAb
Accession(Primary)	Q9Z1B3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 138 kDa; Observed MW: 138 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-:1:200-1:2000

Target information (Q9Z1B3)

Gene Name Plcb1

About Gene Symbol



