

IGF2BP2 Rabbit mAb


ABZ1711

Product Overview

Name	IGF2BP2 Rabbit mAb
Accession(Primary)	Q9Y6M1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66,62 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200


Target information (Q9Y6M1)

Gene Name	IGF2BP2
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of IGF2BP2 in 293T lysates using IGF2BP2 antibody.



Western blot analysis of IGF2BP2 in (1) SH-SY5Y lysates; (2) mouse heart lysates using IGF2BP2 antibody.