

Dnmt3a Rabbit mAb

ABZ1713


Product Overview

Name	Dnmt3a Rabbit mAb
Accession(Primary)	Q9Y6K1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q9Y6K1)

Gene Name DNMT3A

About Gene Symbol



Western blot analysis of Dnmt3a in HeLa, A549 lysates using Dnmt3a antibody.

