


Epsin 1 Rabbit mAb

ABZ1714


Product Overview

Name	Epsin 1 Rabbit mAb
Accession(Primary)	Q9Y6I3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100

Target information (Q9Y6I3)

Gene Name EPN1**About Gene Symbol**

Western blot analysis of Epsin 1 in PC-3, mouse brain lysates using Epsin 1 antibody.



Western blot analysis of Epsin 1 in rat Brain
lysates using Epsin 1 antibody.