

CLIC4 Rabbit mAb


ABZ1716

Product Overview


Name	CLIC4 Rabbit mAb
Accession(Primary)	Q9Y696
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9Y696)


Gene Name	CLIC4
About Gene Symbol	CLIC4 (Chloride Intracellular Channel 4) is a membrane transporter involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CLIC4 in rat Brain, C6, 3T3, Hela lysates using CLIC4 antibody.



Immunocytochemistry analysis of CLIC4(green) in Hela using CLIC4 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CLIC4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.