


RNF115 Rabbit mAb

ABZ1723


Product Overview

Name	RNF115 Rabbit mAb
Accession(Primary)	Q9Y4L5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9Y4L5)

Gene Name RNF115**About Gene Symbol**

Western blot analysis of ZNF364/RNF115 in PC-12 lysates using ZNF364 antibody.



Western blot analysis of ZNF364/RFN115 in Jurkat, HeLa lysates using ZNF364/RFN115 antibody.