

ORP150 Rabbit mAb


ABZ1724

Product Overview

Name	ORP150 Rabbit mAb
Accession(Primary)	Q9Y4L1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 111 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:10-1:100

Target information (Q9Y4L1)

Gene Name	HYOU1
About Gene Symbol	



Western blot analysis of ORP150 in 293T, MCF-7
lysates using ORP150 antibody.