

ORP150 Rabbit mAb


ABZ1725

Product Overview


Name	ORP150 Rabbit mAb
Accession(Primary)	Q9Y4L1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 111 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:50-1:200; ICC/IF-1:50-1:200; FC-1:20; IP-1:20-1:50

Target information (Q9Y4L1)


Gene Name	HYOU1
About Gene Symbol	



Western blot analysis of ORP150 in HeLa lysates using ORP150 antibody.



Immunocytochemistry analysis of ORP150(green) in HeLa using ORP150 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ORP150 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.