


TRAF6 Rabbit mAb

ABZ1727

Product Overview

Name	TRAF6 Rabbit mAb
Accession(Primary)	Q9Y4K3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50

Target information (Q9Y4K3)

Gene Name TRAF6**About Gene Symbol**

Western blot analysis of TRAF6 in 293 lysates
using TRAF6 antibody.

