

TTC11/FIS1 Protein Rabbit mAb


ABZ1729

Product Overview


Name	TTC11/FIS1 Protein Rabbit mAb
Accession(Primary)	Q9Y3D6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100

Target information (Q9Y3D6)


Gene Name	FIS1
About Gene Symbol	



Western blot analysis of TTC11/FIS1 in Jurkat, HeLa lysates using Mitochondrial Fission 1 Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TTC11 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TTC11/FIS1 (green) in MCF-7 using TTC11/FIS1 antibody, and DAPI(blue)