

CAB39 Rabbit mAb


ABZ1733

Product Overview


Name	CAB39 Rabbit mAb
Accession(Primary)	Q9Y376
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (Q9Y376)


Gene Name	CAB39
About Gene Symbol	CAB39 (Calcium Binding Protein 39) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CAB39/MO25 in mouse kidney, rat lung lysates using CAB39 antibody.



Western blot analysis of CAB39/MO25 in Hela, A549, HL-60, U2OS, U87-MG lysates using CAB39/MO25 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CAB39/MO25 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.