

## NDUFAF1 Rabbit mAb


ABZ1734

### Product Overview

|                             |                                                                            |
|-----------------------------|----------------------------------------------------------------------------|
| <b>Name</b>                 | NDUFAF1 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | Q9Y375                                                                     |
| <b>Host Species</b>         | Rabbit                                                                     |
| <b>Clonality</b>            | Monoclonal                                                                 |
| <b>Isotype</b>              | IgG                                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                          |
| <b>Calculated MW</b>        | Calculated MW: 38 kDa; Observed MW: 35 kDa                                 |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:500; FC-1:100-1:200 |

### Target information (Q9Y375)

|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | NDUFAF1 |
| <b>About Gene Symbol</b> |         |



Western blot analysis of NDUFAF1 in rat Brain, MCF-7, C6, 3T3, HeLa lysates using NDUFAF1 antibody.