

NDUFAF1 Rabbit mAb


ABZ1735

Product Overview



Name	NDUFAF1 Rabbit mAb
Accession(Primary)	Q9Y375
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:1000

Target information (Q9Y375)

Gene Name	NDUFAF1
About Gene Symbol	



Immunocytochemistry analysis of NDUFAF1 (green) in LNCap using NDUFAF1 antibody, and DAPI(blue).



Western blot analysis of NDUFAF1 in K562, HeLa lysates using NDUFAF1 antibody.